

SAFETY DATA SHEET

(According to EEC Directive 91/155/EEC and amendments)

ELECTROLYTE SOLUTION FOR K-ELECTRODE F2311K (Code 943-714)

1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Product name: ELECTROLYTE SOLUTION FOR K-ELECTRODE F2311K (Code 943-714)

Application: For in vitro diagnostic use Container size: 15.8 g

Supplier: Radiometer America Inc.
810 Sharon Drive
Westlake OH 44145
440-871-8900 or 800-535-5053 for Emergency

2. HAZARDS IDENTIFICATION

The hazardous properties of the product are considered to be limited.

3. COMPOSITION/INFORMATION ON INGREDIENTS

<u>Chemical name:</u>	<u>CAS-number:</u>	<u>%:</u>	<u>Hazard classification:</u>
Organic and inorganic salts	-	< 2	-

4. FIRST-AID MEASURES

Inhalation: Move into fresh air and keep at rest.

Skin contact: Remove contaminated clothes and rinse skin thoroughly with water.

Eye contact: Flush with plenty of water for up to 15 minutes. Remove any contact lenses and open eyelids widely. If irritation persists: Seek medical attention and bring along these instructions.

Ingestion: Rinse mouth and drink plenty of liquid.

5. FIRE-FIGHTING MEASURES

Extinguishing media: Fire precautions: None specific. Use fire-extinguishing media appropriate for surrounding materials.

Specific hazards: The product is non-combustible.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions: Avoid prolonged contact with skin.

Environmental precautions: Not relevant considering the small amounts used.

Methods for cleaning up: Flush away small quantities with copious amounts of water.
For waste disposal see section 13.

7. HANDLING AND STORAGE

Handling: Observe good industrial hygiene practices.

Storage: Store in closed original container.

SAFETY DATA SHEET

(According to EEC Directive 91/155/EEC and amendments)

ELECTROLYTE SOLUTION FOR K-ELECTRODE F2311K (Code 943-714)

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

<u>Engineering measures:</u>	None specific.
<u>Occupational exposure limits:</u>	Not relevant.
<u>Personal protection:</u>	Protective gloves are recommended.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance:	Liquid
Odour:	Odourless
pH:	5.50
Melting point:	Not available
Boiling point:	Not available
Decomposition temperature:	Not available
Flash point:	Not available
Explosion properties:	Not available
Vapour pressure:	Not available
Vapour density:	Not available
Density:	Abt. 1.0 g/ml
Solubility:	Complete
Other:	Not available

10. STABILITY AND REACTIVITY

<u>Stability:</u>	Stable.
<u>Conditions/materials to avoid:</u>	Not known.
<u>Hazardous decomposition products:</u>	Not known.

11. TOXICOLOGICAL INFORMATION

<u>Inhalation:</u>	Not relevant at normal room temperatures.
<u>Skin contact:</u>	Prolonged contact may cause redness and irritation.
<u>Eye contact:</u>	Splashes may irritate.
<u>Ingestion:</u>	May irritate and cause malaise.
<u>IARC:</u>	Contains small amount of: Phenylmercury nitrate, CAS no. 55-68-5 (IARC Group 3: not classifiable as to carcinogenicity to humans).

12. ECOLOGICAL INFORMATION

<u>Mobility:</u>	Water soluble. Will spread in water systems.
<u>Degradability:</u>	Not known.
<u>Ecotoxicity:</u>	Not known.

SAFETY DATA SHEET
(According to EEC Directive 91/155/EEC and amendments)

ELECTROLYTE SOLUTION FOR K-ELECTRODE F2311K (Code 943-714)

13. DISPOSAL CONSIDERATIONS

EWC-code (European Waste Catalogue): 16 05 06
Abfallschlüsselnummern in Austria (ÖNORM S2100): 59305
Dispose of in accordance with Local Authority requirements.

14. TRANSPORT INFORMATION

Not classified as dangerous goods.

15. REGULATORY INFORMATION

Labelling: An assessment by Radiometer Medical ApS on June 1, 2007 has shown that according to EEC Directives 67/548/EEC and 99/45/EEC, the product is not liable to classification and labelling.

16. OTHER INFORMATION

Limitations/Restrictions: The user must be instructed in the proper work procedure and be familiar with the contents of these instructions.
This SDS complies with the US ANSI standard (ANSI 2400.1-1993) and the international standard (ISO 11014-1)

The information on this data sheet represents our current data and is accurate to the best of our knowledge as to the proper handling of this product under normal conditions and in accordance with the application specified on the packaging and/or technical guidance literature. Any other use of the product which involves using the product in combination with any other product or any other process is the responsibility of the user.

Based on:
Drawing No.: 17479-A4
Version No.: 7

Prepared by:
RADIOMETER MEDICAL ApS

Version No.: 2
Revision date: June 1, 2007