

Glucose for the ABL80 FLEX



Lab-quality glucose in a POC analyzer

Now you don't have to sacrifice quality glucose measurements taken at the point of care. Radiometer's ABL™80 FLEX delivers superior accuracy, stability and interference-free glucose performance in a portable POC analyzer.

Accuracy

- Measuring range of 0-75 mmol/L (1,351 mg/dL)
- Exceptional linearity; largely insensitive to pO_2 levels
- Hematocrit and plasma protein-independent performance

Stability

- Fast start-up
- Solid state – no additional handling or preparation required
- Storage at room temperature

Superior interference properties

- Interference-free performance ($< |0.1|$ mmol/L)
- Demonstrated non-interference from common therapeutic substances

For more information, go to www.radiometeramerica.com/abl80.

Data subject to change without notice.

Radiometer™ and ABL™ are trademarks of Radiometer Medical ApS, Denmark.

RADIOMETER
COPENHAGEN 